

HEK293-6E Cell Preparation for Transfection (20101121JV)

Purpose: The purpose of the guidelines below is to provide a procedure to generate HEK293-6E cells that allow optimal transient protein expression.

1. Materials, Conditions:

Culture media: FreeStyle™ 293 Expression Medium (Invitrogen Cat. 12338-018)

Incubator: 8% CO₂, 37°C, humidified atmosphere. Orbital shaking 90 rpm for 1L culture, 130 rpm for smaller volumes.

Labware: Polycarbonate shake flasks with vented lids. The flasks can be re-used after washing and autoclaving, but have to be checked carefully for potential cracks to avoid media spills in the incubator. The following sizes are used for specific volumes:

Flask Volume (ml)	Culture Volume (ml)	rpm	Corning Cat.#
2800	1000	90	431252
1000	400	130	431147
250	100	130	431144
125	30	130	431143

2. Cell Numbers

The following cell numbers guarantee optimal protein expression

- A. Grow cells at least until they are 2×10^6 but no denser than 3×10^6 .
- B. When splitting the cells, dilute them to no less than 1×10^6 .
- C. Cells need to double with a 24 hr doubling rate to be ready for transfection
- D. Ideally have the cells double two consecutive times at maximum doubling rate (“log phase”) before Transfection
- E. At transfections dilute cells to 1×10^6 .