

Anti FcRn ELISA Protocol

- 1) Coat plate with IG3 diluted 1/1000 (use the 12/1/95 Ab) into PBS. 100 μ l/well, Incubate overnight at 4 C.
- 2) Wash three times with PBS/0.05% Tween pH 7-7.6 (Take excess out by taping the plate against a paper towel).
- 3) Incubate supernatants, 200 μ l/well for 1 hour @ 37 C.
- 4) Wash three times with PBS/0.05% Tween pH 7.0-7.6.
- 5) Make PBS/3%BSA ~pH 7.5

Dilute anti- 2m (polyclonal rabbit a human 2m) into this at 1/1000 dilution. Plate 100 μ l and incubate 1 hour @37 C.

- 6) Wash three times with PBS/0.05% Tween pH 7.0-7.6.
- 7) Dilute peroxidase conjugate -rabbit 1/5000 into PBS/3% BSA, plate 50 μ l/well. Incubate plates 1 hour @ 37 C.
- 7) Wash three times with PBS/0.05% Tween pH 7.0-7.6.
- 8) Prepare a solution containing:

10 mg OPD (ortho phenylene diamine) in
*25 mL of citrate/phosphate buffer pH 5.5
add 12.5 μ l of 30% H₂O₂.

* To prepare Citrate/Phosphate buffer:

20 mL 1M dibasic Sodium Phosphate (pH to 5.5 w/ 1M citric Acid)
170mL distilled water
200mL Final Volume

- 9) Add 200 μ l/well to ELISA plates
- 10) Watch for color development (yellow-->brown, when the solution is yellowish/brownish)
 - Quench rapidly with 50 μ l of 3M sulfuric acid.
 - Read absorbance @ 490 nm.

