

HFE PURIFICATION

Developed by Jose Lebron

Column Info:

Antibody: anti-HFE (heavy chain) mAb 1C3.1G4
Bed volume: 3 - 4 ml
Wild-type HFE capacity: approx. 3 mg
Flow rate: 0.75 ml/min max.
Made by: Jose Lebron (3/96; for info see his notebook #6, p. 176)

Buffers:

<u>A. Binding:</u>	50 mM Tris, pH 7.4 0.05% NaAzide	<u>For 1 L:</u> 50 ml of 1 M stock 5 ml of 10% stock 950 ml dH ₂ O
<u>B. Elution:</u>	50 mM DEA, pH 11 (DiEthylAmine)	<u>For 250 ml:</u> 1.29 ml DEA stock pH with concHCl to 10.9-11 (will create some vapor)

To neutralize fractions: 1 M monobasic NaPhosphate

Column run:

1. Run harvest over column
2. Wash with Buffer A for about 1 hr (about 40 ml)
3. Elute with Buffer B: collect 20 1 ml fractions (25 to elute tail completely). Add 100 ul 1 M monobasic NaPO₄ per tube before elution to neutralize pH. Peak should begin to come off in fraction 8-11, depending on amount of liquid over column.
4. Wash with 20-30 ml Buffer A.
5. Eluted fractions: as long as pH is close to 7, protein may be stored in the DEA/phosphate buffer, but add NaCl to 0.1 M final.

Christina Milburn, 8/26/97

Xtalpower HD\Lab Gen. Info\Protocols\Protein Purif. Protocols\HFE Purification