

PSMA Protein Purification Protocol

Mindy Davis 4-20-05

(1) Protein is obtained from Protein Expression Facility at Caltech. Usually only 1-2 L are ordered since the PSMA expression level is high (>10 mg/L). The supernatants are stored at 4°C.

(2) Concentrate and Exchange

- a) prepare 10kDa cut-off membrane with Amicon 2L concentrator, washing it and pre-equilibrating with 500 mL buffer (50 mM Tris pH=7.5, 150 mM NaCl)
- b) add 0.2 or 0.45 μ m-filtered protein to reservoir, follow concentrator/membrane instructions keeping pressure below 20 psi
- c) concentrate protein to ~200 mL at RT
- d) add 1 L buffer (50 mM Tris pH=7.5, 150 mM NaCl)
- e) repeat steps b)+c) for a total of 3 times (this can also be done by diafiltration for total of 3 L of buffer)
- f) concentrate to 200 mL
- g) drain concentrator into 1L bottle
- h) wash concentrator with ~500 mL buffer to get all protein out into 1L bottle
- i) adjust the ~700 mL of protein so that it is 50 mM Tris pH=7.5, 300 mM NaCl, 10% glycerol, 10mM imidazole (add azide if desired)
- j) mix well
- k) 0.2 μ m filter the protein solution
- l) wash concentrator for next user

(3) Nickel Column

- a) prepare Ni column per Qiagen instructions
- b) use minimal Ni-NTA beads (a mL of beads can bind several mgs of protein)
- c) load O/N at 4°C at rate less than or equal to 1 mL/min
- d) hook column to U/V detector at RT
- e) wash column until baseline is flat at rate less than or equal to 1 mL/min (wash buffer = 50 mM Tris pH=7.5, 300 mM NaCl, 10 mM imidazole, 10% glycerol, azide if desired)
- f) elute column at 2 mL/min (elution buffer = 50 mM Tris pH=7.5, 300 mM NaCl, 250 mM imidazole, 10% glycerol, azide if desired)

(4) SEC

- a) concentrate protein at ~10 °C using amicon concentrator (15 mL size, 10kD membrane)
 - I try to keep protein at less than 8 mg/mL so I will concentrate to 2-4 mL depending on how much supernatant I started with, expression level, etc
- b) pre-equilibrate column (S200 16/60) with wash buffer (50 mM Tris pH=7.5, 150 mM NaCl)
- c) I do 2 mL or smaller injections onto this column so sometimes I will do several runs if I have a lot of protein
- d) run column at 1 mL/min and collect fractions (1.5 mL/fraction) using Akta program
- e) run gel of fractions containing protein (7.5 % SDS-Page Reducing conditions; PSMA runs near 97.4 MW marker)
- f) the dimeric PSMA fractions are what I use for crystallography applications

- (5) If needed: Factor Xa cleavage will remove the N-terminal 6x-His-tag
- a) if ~10 mL of protein from the sizing column use 2.5 mL water, 20 μL CaCl_2 (from 1M stock) and 40 μL of Factor Xa. (this protocol was from another protein in the lab and I didn't optimize it; the Factor Xa instructions for their assay conditions are 1 μg Factor Xa to 50 μg of their test protein. Their reaction was in 50 μL , 20 mM Tris-HCl, 100 mM NaCl, 2 mM CaCl_2 (pH=8.0))
 - b) this reaction is allowed to go for 24 hours at 4°C on a very slowly turning rotator (very fast rotation leads to extra ppt)
 - c) the protein is purified via Ni-NTA (removes any uncleaved or free tags) and the SEC as above to separate Factor Xa from PSMA. (There is a Factor Xa you can buy that binds to Ni-NTA but I haven't tried that)