

## Gene Construction (via Restriction Cloning):

**Step 1:** Design and order oligos (25nmol scale, desalted)

**Step 2:** Resuspend each lyophilized oligo to 1  $\mu\text{g}/\mu\text{L}$  with sterile water

**Step 3:** Make 40  $\mu\text{L}$  of a 10  $\mu\text{M}$  stock of each oligo

**Step 4:** Make the PCR reaction mix:

5  $\mu\text{L}$  plasmid template (10 ng/ $\mu\text{L}$ ;  $\sim 50$  ng)  
1  $\mu\text{L}$  fwd primer (10  $\mu\text{M}$ )  
1  $\mu\text{L}$  reverse primer (10  $\mu\text{M}$ )  
18  $\mu\text{L}$  water  
25  $\mu\text{L}$  Q5 Polymerase 2X MM (NEB; or other MM)

If using a standard polymerase, make sure to add 10X buffer and dNTPs.

**Step 5:** Run PCR reaction (72°C extension temperature for Q5)

95°C	2 min		
95°C	30 sec	}	10 cycles
55°C	30 sec		
72°C	30 sec		
95°C	30 sec	}	20 cycles
60°C	30 sec		
72°C	30 sec		
72°C	10 min		
hold at 4°C			

**Step 6:** Run the PCR Cleanup Kit (Qiagen) on the PCR product

- Elute the DNA from the column with 50  $\mu\text{L}$  of water

**Step 7:** Set up the restriction digestion of the PCR product and the vector

- Consult the NEB (New England Biolabs) catalog or website for the appropriate reaction conditions (buffer and BSA)
- Set up 20 $\mu\text{L}$  reactions (each enzyme individually—controls, vector dd and PCR dd)
  - Pipette each reaction gently to resuspend the enzyme and mix completely
- Incubate at 37°C for 2 hours

### Enzyme 1 only

2  $\mu\text{L}$  10X Buffer  
X  $\mu\text{L}$  vector plasmid ( $\sim 1000$  ng)  
X  $\mu\text{L}$  water (20 $\mu\text{L}$  total rxn)  
1  $\mu\text{L}$  enzyme 1

### Enzyme 2 only

2  $\mu\text{L}$  10X Buffer  
X  $\mu\text{L}$  vector plasmid ( $\sim 1000$  ng)  
X  $\mu\text{L}$  water (20 $\mu\text{L}$  total rxn)  
1  $\mu\text{L}$  enzyme 2

### Vector dd

2  $\mu\text{L}$  10X Buffer  
X  $\mu\text{L}$  vector plasmid ( $\sim 2000$  ng)  
X  $\mu\text{L}$  water (20 $\mu\text{L}$  total rxn)  
1  $\mu\text{L}$  enzyme 1  
1  $\mu\text{L}$  enzyme 2

## PCR Product dd

2  $\mu$ L 10X Buffer  
14  $\mu$ L Amp PCR product  
1  $\mu$ L enzyme 1  
1  $\mu$ L enzyme 2

### Optional Step 7b: Dephosphorylate vector

- Calf Intestinal Alkaline Phosphatase (CIP) is active in CutSmart Buffer
- Add 1 $\mu$ L CIP to your VECTOR dd reactions -- do not add to your insert dd tubes; mix by pipetting
- Incubate at 37°C for 30-60 minutes

### Step 8: Gel purify double-digested (dd) vector and insert

- Pour a 1% Agarose gel (1X TAE) - 50mL for small gel; 120 mL for large gel
- Add 10X DNA Loading buffer to each of the samples
  - Dark blue (bromophenol blue) runs at ~500 bp  $\rightarrow$  use this for large DNAs (vector, etc)
  - Light blue (xylene cyanol) runs at ~3000 bp  $\rightarrow$  use this for small DNAs (PCR products, etc)
- Load the gel
  - 1  $\mu$ L NEB 2-log ladder + 8  $\mu$ L water + 1  $\mu$ L BPB dye (dark blue)  $\rightarrow$  Load all 10  $\mu$ L
  - ~1000 ng uncut vector
  - Load ALL of each digestion reaction
  - Load uncut PCR product if desired
- Run for 45-70 min at 100 V (for TAE buffer) until the bromophenol blue runs ~3/4 of the way down the gel
- Extract dd vector and insert from the gel and run the Qiagen Gel Extraction kit on each band (elute from column in water)

### Step 9: Ligate the dd PCR product to the dd vector

- Determine the concentration of double-digested vector and insert from the gel purification by nanodrop
  - Use single digestion controls to determine the % of cleavage
- Determine the amount of insert that is needed for a 1:3 molar ratio (vector: insert)
  - This equation will tell you how much insert to use for a 1:1 molar ratio, multiply by 3 to get the amount for the 1:3 ligation

$$\frac{Xbpvector}{Xbpinsert} = \frac{50ngvector}{Xnginsert}$$

#### Control Ligation

1  $\mu$ L 10X T4 Ligase Buffer  
X  $\mu$ L dd vector (~50 ng)  
X  $\mu$ L water  
1  $\mu$ L T4 DNA Ligase

#### 1:3 Ligation

1  $\mu$ L 10X T4 Ligase Buffer  
X  $\mu$ L dd vector (~50 ng)  
X  $\mu$ L dd insert (~1:3 molar ratio)  
X  $\mu$ L water  
1  $\mu$ L T4 DNA Ligase

- Pipette the reactions gently to mix
- Incubate ligation reactions at room temperature for 1-2 hours

### Step 10: Transform ligation reactions into Z-comp cells

- Thaw appropriate molecular biology Z-comp cells (XL1-blue, TOP10, etc) on ice (1 tube contains 100  $\mu$ L of cells—enough for 2 transformations)
- Add 2  $\mu$ L of each ligation into 50  $\mu$ L of cells, stir with pipette tip (do not pipette up and down or vortex), let sit on ice for 2-3 minutes
- Plate entire transformation onto a warmed agar plate
- Incubate overnight at 37°C

### Step 11: Count colonies on plates, record number of colonies on vector-insert ligation plates vs. control ligation plates

### Analyze your results if possible before sequencing: Colony PCR or Restriction Analysis

#### 1. Colony PCR:

- **Day 1:** Determine how many colonies to do PCR on (If there are many more colonies on the ligation than the control ligation, 6-10 colonies is sufficient) – determine based on how many lanes you can run in a gel
- Circle and number the colonies on the plate
- Control(s):
  - Colony from control ligation plate
  - Plasmid that was used for vector (~5-10 ng of purified plasmid/rxn)
  - Positive control plasmid if available
- Prepare Master Mix of reaction (30  $\mu$ L/reaction – I usually make one more volume than necessary):

#### Colony PCR

3  $\mu$ L 10X Thermopol Buffer  
0.3  $\mu$ L 10mM dNTPs  
1  $\mu$ L fwd primer (10  $\mu$ M)  
1  $\mu$ L rev primer (10  $\mu$ M)  
24.4  $\mu$ L water  
0.3  $\mu$ L Taq polymerase

- Distribute 30  $\mu$ L of reaction mixture into each PCR tube
- With a sterile toothpick carefully pick a colony and place the toothpick in the PCR tube
- Add plasmid DNA to any controls where it is needed
- Run PCR reaction

95°C	2 min	} 25 cycles
95°C	30 sec	
50°C	30 sec	
68°C	1 min/kb	
68°C	5 min	
hold at 4°C		

- Add gel running buffer (3  $\mu$ L/samples) and run the completed reactions on a 1% Agarose gel (load all 30  $\mu$ L)
- Pick the remaining colony into 3mL LB + Amp and grow overnight for minipreps
- **Day 2:** Miniprep and send for sequencing

## 2. Restriction Analysis:

- **Day 1:** Determine how many colonies to do restriction analysis on (If there are many more colonies on the ligation than the control ligation, 6-10 colonies is sufficient) – determine based on how many lanes you can run in a gel
- Pick colonies into 3mL LB + Amp and grow overnight for minipreps
- **Day 2:** Miniprep cultures
- Set up restriction analysis (10uL reactions)
  - Choose an enzyme or a pair of enzymes that give you a different product for your desired construct compared to the template you started with.
  - Examples:
    - Use an enzyme that cuts once in your template but twice in your desired construct; look for a small digested fragment
    - Use the same set of enzymes you used to do your cloning if you are inserting a different sized fragment than you started with (linker cloning)
  - Make sure to run appropriate controls:
    - Single cut on template construct
    - DD on template construct
    - DD on colony from control ligation plate
  - Make a master mix of the dd and add to 1-2 uL of miniprep DNA (make 1 more MM volume than you expect to use)

### **dd MM**

1  $\mu$ L 10X Buffer  
6  $\mu$ L water  
0.5  $\mu$ L enzyme 1  
0.5 uL enzyme 2

- Incubate restriction digestions for 30-60 minutes
- Run a 1% agarose gel
  - Make sure to run uncut template vector, single cut controls, and control colony controls.
- Select positive constructs and send to sequencing