

Sample Preparation for Mass Spectroscopy using ZipTip_{C18}

1. Make solutions.

Sample preparation solution (2.5% trifluoroacetic acid)

Combine: 25 μ l trifluoroacetic acid
475 μ l nanopure water

Wetting solution (50% acetonitrile)

Combine: 500 μ l acetonitrile
500 μ l nanopure water

Equilibration solution (0.1% trifluoroacetic acid)

Combine: 1 μ l trifluoroacetic acid
999 μ l nanopure water

Wash solution (0.1% trifluoroacetic acid)

Combine: 1 μ l trifluoroacetic acid
999 μ l nanopure water

Elution solution (50% acetonitrile, 0.1% trifluoroacetic acid)

Combine: 500 μ l acetonitrile
499 μ l nanopure water
1 μ l trifluoroacetic acid

2. Dilute the purified protein in water*, **final amount: 20 μ l of a 50 μ M solution.** This is a higher concentration than the facility requests, however I have observed that our proteins do not bind well to the ZipTip matrix. Therefore, I overload the tip to compensate and that has produced good Mass Spec data.

*Maximum binding to the ZipTip_{C18} is achieved in the presence of TFA or other ion-pairing agents. To maximize analyte binding, the **final TFA concentration should be between 0.1% and 1% at a pH of <4.0.**

3. Prewet the ZipTip_{C18} by slowly aspirating 10 μ l of Wetting solution into the ZipTip_{C18}. Dispense to waste. Repeat.
4. Equilibrate the tip for binding by washing it twice with 10 μ l of Equilibration solution.
5. Bind the protein to the ZipTip_{C18} by slowly aspirating and dispensing the sample 10 times.
6. Wash the bound protein three times with Wash solution by slowly aspirating and dispensing to waste.
7. Dispense 10 μ l of Elution solution into a clean tube using a standard pipet tip.
8. Carefully, aspirate and dispense the 10 μ l of Elution solution through the ZipTip_{C18} five times, **without introducing air into the sample.**
9. Submit for Mass Spec!