

Restriction of PCR product:

If you had a bright DNA band after PCR you can do a cut on 20 ul PCR product. If you had a faint band, speed vac the 50 ul you got from the gel extraction and use 20 ul concentrated product.

- do the restriction in a total 50 ul mix
- choose the buffer that is good for both enzymes you are going to use and add 5 ul
- add 1 ul of each enzyme
- incubate at 37° C (or other temperature) for 1 hour or longer
- heat inactivate 20' at 65° C (if the enzymes are heat inactivatable)

Restriction of vector:

- cut 5-10 ug of vector in a 50 ul mix for 1.5 hours using the first enzyme
- run the DNA on a 1% LMP agarose gel
- cut out the (cut) DNA band and gel extract
- take 25 ul and cut this with the second enzyme in a 50 ul mix for 1.5 hours save the other 25 ul
- run the DNA on a 1% LMP agarose gel again and cut out the DNA band

Ligation

* mix the following things together:

- 1 ul cut vector
- 5 ul insert
- 1 ul Mg⁺⁺ (20x)
- 5 ul H₂O

do a negative control (only vector)

- mix
- incubate 5' at 37° C
- add 2 ul 10x ATP/DTT (10mM/100mM)
- add 1 ul T4 DNA ligase
- incubate 2 hours at room temperature then ON at 16° C

Transformation:

Use electrocompetent cells with an efficiency of at least 10e7

- thaw the competent cells on ice
- add 1 ul ligation mix
- electroporate at 25 uF, 200 Ω , 2.5 Kvolt
- add 1 ml of SOC media and incubate 1 hour in the shaker at 37° C
- spin cells down 6 sec at 14000 rpm
- take ~ 800 ul SOC out and resuspend the bacteria in the left over media
- plate cells on a agarose plate containing the right antibiotic
- let grow ON in 37° C incubator