

Favored *E. coli* strains for site-directed mutagenesis using uracil replacement:

RZ1032: dut<sup>ts</sup>, ung, tet<sup>r</sup>

Bw313:

**dut:** deoxy uracil transphosphatase (enzyme required for formation of thymidine precursor)

**ung:** uracil-DNA-glycosylase (enzyme which normally will proof and replace uracil in DNA)

## STEP 1. GEL PURIFY MUTAGENESIS OLIGOS (IF REQUIRED)

- 1a. Purify crude oligos from 1% agarose gel.
- 1b. Desalt oligo prep using NAP 25 column.
- 1c. EtOH precipitate oligo DNA and resuspend into TE to give final concentration of 1.8 µg/ .

**NOTE:** For ss DNA, 1 OD<sub>260</sub> = 20 µg/ml.

## STEP 2. ALTERNATIVE OLIGO KINASE REACTION (20λ total volume)

- 2a. Add the following:

13	dH <sub>2</sub> O
2	oligo DNA (approx. 10 ng)
2	10X T4 ligase buffer (Boeh. Mann.)
1	10mM ATP
<u>2</u>	T4 kinase (Boeh. Mann.)
20	

- 2b. Incubate for 3 hr at 37°C.
- 2c. Store kinased oligos at -70°C until ready for use.

## STEP 3. PREPARATION OF TEMPLATE SINGLE STRAND DNA (ss DNA)

- 3c. Prepare Ca<sup>2+</sup> competent RZ1032 cells (these cells can't be made electrocompetent).
- 3d. Transfect template (ds DNA) into Ca<sup>2+</sup> competent RZ1032 cells:
  - add 0.1 µg of template ds DNA to a cold eppendorf tube.
  - add 50 of thawed CA<sup>2+</sup> comp. RZ1032 cells, invert to mix reaction contents.
  - incubate reaction 5 min on ice, then 5 min at 37°C.
  - immediately add 1 ml of 2XYT media and incubate reaction 1 hr at 37°C.
  - streak reaction Hi/Lo on LB + Amp + Tet plates, place plate in 37°C incubator until colonies appear.
- 3e. Pick a single colony from transfect reaction plate and grow in 2 ml of LB (or 2XYT) media containing 100 µg/ ml Amp, 10 µg/ ml Tet, and 0.25 µg/ml fresh Uridine in a 37°C shaker incubator.

- 3f. After approximately 8 hrs (usually overnight, O/N), removed 100 of the shaker culture and add 2 ml of 2XYT media containing 20 of helper phage (either M13K07 or VCSM13 strain), 10 µg/ml Tet and 0.25 µg/ml fresh Uridine (no Amp or Tet is added to media at this time).  
**NOTE:** assume 1 µg of ss DNA from 1 ml of phage culture, e.g., to recover 10µg of ss DNA start with 10 ml of phage culture (5x 2 ml cultures). Remember, one site-directed mutagenesis reaction requires a minimum of 5µg of ss template DNA. If using phage culture larger volumes (10-20 ml), use 50 ml SS-34 rotor tubes with yellow caps for incubating (see step 3j).
- 3g. Incubate phage culture 1 hr in 37°C shaker incubator.
- 3h. After 1 hr incubation, add to phage culture Amp (final conc. = 50 µg/ml) and Kanamycin (final conc. = 50 µg/ml).
- 3i. Incubate phage-antibiotic culture for 8 hrs in 37°C shaker incubator.
- 3j. Pellet the phage-antibiotic culture at 1900xg for 10 min (4000 rpm using SS-34 rotor).
- 3k. Pool supernatants from multiple culture tubes into one sample. Add 1 vol. of 20% PEG, 2.5 M NaCl for each 4 vol. of supernatant. Mix by inverting tube and incubate reaction overnight at 4°C.
- 3l. Pellet PEG/NaCl reaction at 7650xg for 20 min (8000 rpm using a SS-34 rotor).
- 3m. Remove and discard supernatant. The pellet containing the precipitated phage particles should be clearly present at base of the rotor tube. In addition, there may a slight whitish film of PEG lining the rotor tube wall. Using a very small amount of dH<sub>2</sub>O (200 ), wash the rotor wall covered by the PEG film until the PEG dissolves completely into dH<sub>2</sub>O wash. Repeated washings may be required to remove the PEG film. Be careful to not disturb the phage pellet.
- 3n. Resuspend the phage pellet with 500 TE pH 8.0 and transfer suspension to 1.5 ml eppendorf tube. Do not use a TE with a lower pH as your ss DNA will dissolve into the phenol phase in the extraction procedure following this step.
- 3o. Extract the phage pellet resuspension with an equal volume of buffer phenol pH 8.0. Remove the aqueous (top) phase containing the ss DNA and repeat 1:1 volume phenol extraction until the interface does not contain a whitish ppt. Follow the phenol extractions with one chloroform extraction. Add two volumes of chloroform for each vol. of ss DNA in TE.
- 3p. Precipitate the ss DNA by adding 1/10 vol. of 3M NaAc and 2 vol. of 100% EtOH. Place reaction in -70°C freezer for a minimum of 2 hrs.
- 3q. Pellet the precipitate ss DNA at 12000 rpm, 4°C for 30 min using microfuge.
- 3r. Wash ss DNA pellet using 80% EtOH and dry the pellet by speedvac.  
**NOTE:** place parafilm across 1.5 ml tube opening to prevent pellet being lost when releasing speedvac vacuum.
- 3s. Resuspend the ss DNA pellet in 50 TE.
- 3t. Measure concentration and purity of ss DNA by OD<sub>260/280</sub> ratio. Verify quality and quantity of ss DNA by running a sample on a 1% agarose gel containing a known quantity of ss DNA.

**\*STEP 4. ALTERNATIVE TEMPLATE-PRIMER ANNEALING REACTION (10 $\lambda$  total volume)**

4a. Add the following:

2     autoclaved ddH<sub>2</sub>O  
6     template ss DNA (5  $\mu$ g)  
1     kinased oligo (500ng)  
1    10X T4 DNA ligase buffer (NEB)  
10

4b. Incubate 70°C for 2 min, chill on ice for 2 min, incubate at 22°C (RT) for 10 min.

**\*STEP 5. ALTERNATIVE ELONGATION/LIGATION REACTION (100 $\times\lambda$  total volume)**

5a. Add the following to the 10 alternative annealing reaction:

40    dH<sub>2</sub>O  
25    2mM dNTP  
10    10mM ATP  
10    10X T4 DNA ligase buffer (Boeh. Mann.)  
3     T4 DNA polymerase (Boeh. Mann.)  
2    T4 DNA ligase (Boeh. Mann.)  
90  
  
+ 10    alternative annealing reaction  
100    total volume, elongation/ligation reaction

5b. Incubate reaction at 22°C (RT) for 2 hrs.

5c. Incubate reaction at 37°C for 1 hr.

5e. Check result of elongation/ligation reaction by running a sample on a 1% agarose gel. Include the following controls:

- i) ss DNA, no primer, no elongation
- ii) ss DNA, no primer, but elongated
- iii) ds DNA which is equal expected kilobase length of ds mutagenesis DNA product

## STEP 6. PRECIPITATE AND RECOVER MUTAGENIC DS DNA

- 6a. Precipitate ds DNA from elongation/ligation reaction. For each 100  $\mu$ l of elongation/ligation reaction, add 12.5  $\mu$ l of 7.5 M Ammonium Acetate and 375  $\mu$ l of 100% EtOH. Incubate precipitate reaction at -80°C for a minimum of 1 hr.
- 6b. Pellet precipitated ds DNA at 12000 rpm for 30 min using microfuge.
- 6c. Wash ds DNA pellet using 80% EtOH and dry the pellet by speedvac.  
**NOTE:** place parafilm across 1.5 ml tube opening to prevent pellet being lost when releasing speedvac vacuum.
- 6d. Resuspend the ds DNA pellet in 20  $\mu$ l TE.
- 6e. Electroporate 1  $\mu$ l of ss DNA with 50  $\mu$ l of electrocompetent TGI cells.
- 6f. Streak electroporation reaction Hi/Lo on LB + Amp plates.

**STEP 2. OLIGO KINASE REACTION (20 $\lambda$  total volume)**

\* *SEE ALTERNATIVE PROTOCOL*

2a. Add the following:

- 13 dH<sub>2</sub>O
- 1 oligo DNA (final conc. = 250ng/ )
- 2 10X Salt Stock (500mM Tris pH 7.8, 100mM MgCl<sub>2</sub>)
- 2 10X ATP/DTT (10mM ATP, 100mM DTT)
- 2 T4 kinase (NEB: New England Biolab)
- 20

2b. Incubate for 3 hr at 37°C.

2c. Store kinased oligos at -70°C until ready for use.

**STEP 4. TEMPLATE-PRIMER ANNEALING REACTION (20 $\lambda$  total volume)**

\* *SEE ALTERNATIVE PROTOCOL*

4a. Add the following:

- 6 dH<sub>2</sub>O
- 10 template ss DNA
- 2 kinased oligo (500ng)
- 2 10X Salts (500mM Tris, 100mM MgCl<sub>2</sub>)
- 20

4b. Prepare ss DNA, no primer, but elongated negative control reaction by replacing kinased oligo volume with equal volume of dH<sub>2</sub>O.

4b. Incubate 70°C for 2 min, chill on ice for 2 min, incubate at 22°C (RT) for 10 min.

**STEP 5. ELONGATION/LIGATION REACTION (105 $\lambda$  total volume).**

\* *SEE ALTERNATIVE PROTOCOL*

5a. Add the following to the 20  $\mu$ l annealing reaction:

35	dH <sub>2</sub> O
25	2mM dNTP
10	ATP/DTT stock (10mM ATP, 100mM DTT)
10	10X Salts stock
3	T4 DNA polymerase (Boeh. Mann.)
<u>2</u>	T4 DNA ligase (Boeh. Mann.)
85	

+ 20  $\mu$ l annealing reaction

105  $\mu$ l total volume, elongation/ligation reaction

5b. Incubate reaction at 16°C for 10 min.

5c. Incubate reaction at 22°C (RT) for 2 hrs.

5d. Incubate reaction at 37°C for 1 hr.

5e. Check result of elongation/ligation reaction by 1% agarose gel electrophoresis. Include the following controls:

i) ss DNA, no primer, no elongation

ii) ss DNA, no primer, but elongated

iii) ds DNA which is equal expected kilobase length of ds mutagenesis DNA product