


StrepII-tag Purification

The low affinity of the strepII-StrepTactin interaction ($k_d \sim 1 \mu\text{M}$) causes migration of the adsorption profile down the column (see figure below) in situations of low or zero concentration, as found in dilute samples and during the wash step. *Lengthening the column by connecting 2-3 columns in series and lowering the wash volume prevents the target protein from eluting off the column.*



Ind. Eng. Chem. Res. **2002**, *41*, 2320-2329

Chromatographic parameters	
Column	StrepTrap HP 5mL, 3 in series
Binding Buffer	100mM Tris 150mM NaCl 1mM EDTA pH 8.0
Elution Buffer	100mM Tris 150mM NaCl 1mM EDTA 2.5mM desthiobiotin pH 8.0
Gradient	0-2.5mM desthiobiotin step
Sample Loading Volume	As low as possible (ideally < 100 mL)
Sample Loading Flow Rate	2.5 mL/min
Flow Rate	2.5 mL/min
Wash Volume	1.5 CV (22.5 mL)

1. *Sample preparation.* Concentrate the clarified supernatant 10-20x and buffer exchange into TBS pH 7.4 using tangential flow filtration. The sample volume should be as low as possible. The buffer exchange should be done using at least a 1:10 dilution and then reconcentrated. Residual biotin from the 293 media (Freestyle 293 contains 98 $\mu\text{g/L}$) binds irreversibly to the StrepTactin ligand which will interfere with binding of the target protein and can inactivate the column. (Addition of avidin was found to be unnecessary.)
2. *Column Equilibration.* Equilibrate the column with 5 CV binding buffer, or until baseline is reached in UV and conductivity. Note: All buffers and samples should be 0.45 μM filtered and degassed.
3. *Sample Loading.* Load the sample onto the column at 2.5 mL/min. For large sample volumes, load half of the sample onto the 3 columns in series, then disconnect the 1st column (entrance). Continue loading onto columns 2 and 3 and then reconnect the 1st column. Collect and save the flow through and save a small volume of the load for SDS-PAGE analysis.

4. *Wash.* After sample application, wash out unbound sample with 1.5 CV (22.5 mL) of binding buffer
5. *Elution.* Elute the strepII-tagged protein by competitive elution using a 2.5mM desthiobiotin step gradient for 2 CV (30 mL).
6. *SDS-PAGE analysis.* Run a gel on the load, flow through, and eluate to check for recovery and purity.
7. *Regeneration.* Remove the bound desthiobiotin by using the following regeneration procedure:
 - 3 CV water
 - 3 CV 0.5M NaOH
 - 3 CV water
8. Repeat steps 2-7 on the flow through of the first purification if SDS-PAGE indicates there is still significant protein remaining in the flow through.

The 5 mL columns should be rotated after each use to prevent a rise in backpressure on the first column. The top column should be moved to the bottom (i.e. from 1-2-3 to 2-3-1 to 3-1-2 to 1-2-3).

Checking the integrity of the column:

HABA provides a visual check of column integrity. HABA undergoes a change in color from yellow to red upon complex formation with StrepTactin. If the upper zone does not turn red, this indicates the column has become inactivated.

- 15 CV Binding Buffer + 1mM HABA
- 30 CV Binding Buffer

Cleaning-in-place (CIP) procedure

If there is an increase in backpressure that does not lower during regeneration, the following cleaning procedure can be used:

- 3 CV water
- 3 CV 0.5M NaOH (removes precipitated proteins)
- 3 CV water
- 3 CV 20% EtOH (removes hydrophobically bound proteins)
- 3 CV water

If the nature of the contaminant(s) is known, the appropriate chemical can be used for CIP as long as it is compatible. The compatibility guide can be found at [http://www.gelifesciences.com/aptrix/upp00919.nsf/Content/FDE4EB35510A8A6CC1257628001D5EAD/\\$file/28913630AB.pdf](http://www.gelifesciences.com/aptrix/upp00919.nsf/Content/FDE4EB35510A8A6CC1257628001D5EAD/$file/28913630AB.pdf).

References

Schmidt, T.G., Skerra, A. The Strep-tag system for one-step purification and high-affinity detection or capturing of proteins. *Nat. Protocols* **2**, 1528-1535 (2007)

B. Cass et al., Purification of recombinant proteins from mammalian cell culture using a generic double-affinity chromatography scheme. *Protein Expr. Purif.* **40** 77-85 (2005)