

## Purification of Antibodies

### Protein A:

The buffers I use are the ones recommended by Pharmacia in their catalog. There are several Pharmacia HiTrap protein A columns in the lab; if one of these is used, the column can be run on the FPLC. If you are using a column that has been used for purifying an antibody other than the one you wish to purify, **thoroughly** wash the column with elution buffer to remove any old antibody that may have been left behind. I typically will run 100 ml elution buffer over a column and equilibrate it with 5-10 column volumes of wash buffer before I use it. The wash buffer is 20 mM NaPO<sub>4</sub>, pH 7.0 with 0.05% azide; and the elution buffer is 0.1 M citric acid-NaOH, pH 3.0, also with 0.05% azide. Run your antibody sample over the column, wash with 5-10 column volumes wash buffer, and elute with 5-10 column volumes elution buffer. Collect the eluate into 1 ml fractions; the collection tubes should contain 100  $\mu$ l 1 M Tris, pH 7.4 to neutralize the acid. Concentrate the fractions containing protein (dialyze into an appropriate buffer, e.g. PBS/0.05% azide) and check the concentration by BCA assay and store at 4°C. Wash the column with 5-10 column volumes of wash buffer and run more protein over it. Note: the elution peak from the column starts out sharply, but diminishes very slowly. Pharmacia claims the binding capacity of the HiTrap column is 20 mg IgG/ml of gel. I have never tested the column to see if this is true.

### Protein G:

Exactly the same as protein A purification, except Pharmacia says to use 0.1 M glycine-HCl, pH 2.7 for the elution buffer.

### FcRn column:

The FcRn column must be used in the cold room. The protein sample must first be brought to pH 6.0 with NaPO<sub>4</sub> (bring the concentration of NaPO<sub>4</sub> to 50 mM in the sample). Run the sample over the column; wash with 5-10 column volumes of 50 mM NaPO<sub>4</sub>, pH 6.0; and elute with NaPO<sub>4</sub>, pH 8.0. All the buffers should contain 0.05% azide. Collect the eluate into 1 ml fractions and concentrate those containing protein (dialyze into an appropriate buffer, e.g. PBS/0.05% azide). Check the concentration by BCA assay and store at 4°C.

## Antibodies I Have Purified

### 34-1-2S:

This is a mouse IgG<sub>2a</sub> that recognizes H-2 K<sup>d</sup> and H-2 D<sup>d</sup> MHC class I heavy chains. Purify on protein A columns. Warning: do not give this to Susan Ou to produce ascites without telling her that it will kill mice with the MHC haplotypes: K<sup>d,b,s,r,q,p</sup> and D<sup>d</sup>.

### M1/42.3.9.8.HLK (M1/42 for short):

This is a rat IgG<sub>2a</sub> that reacts with all H-2 haplotypes if they are complexed with mouse  $\beta_2m$ . This MAb cannot be purified on protein A or G columns, and must be purified on the FcRn column.

### H57-597 (H57 for short):

This is a hamster IgG that recognizes the  $\alpha$  subunit of T cell receptors. Purify on a protein G column.

### D4:

This antibody recognizes I-E<sup>k</sup> complexed with a peptide. Purify on a protein A column.