

***E.coli* Expression and Renaturing Protocols**

1. High efficiency competent *E.coli* cells

- (i) Grow bugs in LB or similar - 500 ml for good batch
- (ii) When grown to OD 600 of 0.6, spin down gently.
- (iii) Resuspend in 100mls cold TFB1 and then spin down again
- (iv) Resuspend in 20mls cold TFB2 - aliquot and store at -80°C.

TFB1

30mM KOAc
50mM MnCl₂
100mM KCl
10mM CaCl₂
15% v/v glycerol

TFB2

10mM Na-MOPS pH 7.0
75 mM CaCl₂
10mM KCl
15% v/v glycerol

2. Expression media

	Low salt LB, 1 litre	Normal LB, 1 litre	TYP, 1 litre
Bactotryptone	10g	10g	16g
Yeast extract	5g	5g	16g
NaCl	5g	10g	5g
K ₂ HPO ₄ (anhydrous)	-	-	2.5g

3. Inclusion body preparation

Triton Wash

50 mM Tris-HCl pH 8.0	50mls 1M
0.5% Triton-X100	5mls 100%
200 mM NaCl	40mls 5M
10 mM EDTA	20mls 0.5M
0.1 % (w/v) Sodium Azide	10mls 10%
2 mM DTT, pH 8.0	1ml 2M (or 100ul/100mls)

Make up to 1 litre with water

Resuspension wash

50 mM Tris-HCl pH 8.0	50mls 1M
1 mM EDTA	2mls 0.5M
0.1 % (w/v) Sodium Azide	10mls 10%
2 mM DTT	1ml 2M (or 100ul/100mls)

Make up to 1 litre with water

4. Inclusion body solubilisation

Urea solution

8M Urea	40 mls 10M urea
50mM MES pH 6.5	2.5 mls 1M Mes pH 6.5
10mM EDTA	1 ml 0.5M EDTA
2mM DTT	50ul 2M DTT
H ₂ O	6.45 mls

total: 50 mls

Guanidine solution

6M GdmHCL	15mls 8M Gdm
50mM Mes pH 6.5	1ml of 1M Mes pH 6.5
10mM EDTA	400ml of 0.5M EDTA
2mM DTT	20ml of 2M DTT
H ₂ O	3.58 mls

total: 20mls

4. Refolding buffer

100mM Tris pH 8.0	100mls 1M
2mM EDTA	4mls 0.5M
0.4M L-Arginine-Hydrochloride	84.28g
5mM Oxidised Glutathione	0.31g
0.5mM Reduced Glutathione	1.54g
0.1mM PMSF	1 ml 0.1M

Make up to 1litre with water

5. FPLC buffer

150mM NaCl	30mls 5M
20mM Tris pH 8.0	20mls 1M

Make up to 1 litre with water