

Expression Construct: _____ Cell type _____

Day 1, late PM (date = _____)

1. Grow overnight culture at 37°C, with shaking. Inoculate at _____ PM

Culture volume/container size:

Broth:

Additives:

Day 2 (date = _____)

1. Harvest overnight culture at _____ AM, pellet (1400 x g, 4°C, 10 min) and resuspend in _____ ml fresh LB broth.

2. Inoculate new culture(s) with _____ ml of resuspended cells (dilution = _____), grow at 37°C. Time inoculated: _____ AM
_____ (sample 1) _____ (sample 2) _____ (sample 3) _____ (sample 4)

Culture volume/container size:

Broth:

Additives:

3. Make 100 mM IPTG stock, and keep on ice. (_____ mg in _____ μ l PBS).

4. Remove 1 ml pre-induction sample for SDS-PAGE analysis. Pellet cells (11,000 x g, 5 min, RT), remove SN, resuspend in 100 μ l PBS, add 100 μ l sample buffer – freeze.

5. Induce with IPTG. _____ (sample 1) _____ (sample 2) _____ (sample 3) _____ (sample 4)

OD₆₀₀ at induction:

Time induced:

IPTG volume added: _____ μ l of 100 mM IPTG to each container (giving a final concentration of _____ mM IPTG).

6. Remove 1 ml post-induction sample for SDS-PAGE analysis. Pellet cells (11,000 x g, 5 min, RT), remove SN, resuspend in 100 μ l PBS, add 100 μ l sample buffer – freeze.

OD₆₀₀ at end of induction: _____ (sample 1) _____ (sample 2) _____ (sample 3) _____ (sample 4)
(1:2 with LB)

Time harvested: _____ PM (Total length of induction _____ hrs). Average OD₆₀₀ at harvest: _____

7. Store cultures at 4°C, O/N.

Day 3 (date = _____)

1. Pool and pellet cultures (6,000 x g, 10 min, 4°C). Culture volume harvested: _____. (\div _____ = _____ ml/bottle).

2. Remove, neutralize (10% bleach-final concentration, let sit 1 hour under hood), then discard supernatant.

3. Resuspend cells in ice-cold 20 mM Tris, pH 8.0, 20% sucrose, to an OD₆₀₀ of 5 U/ml, by SWIRLING (no pipeting!). This takes ~20 min.

Volume: _____ ml/pellet

4. Once cells are completely resuspended, incubate on ice (with gentle shaking), 10 min.

5. Pellet cells (15,000 x g, 30 sec, 4°C).

6. Remove, neutralize and discard supernatant.
7. GENTLY resuspend cells by SWIRLING (no pipeting!) in ice-cold 15 mM Tris, pH 8.0, to an OD₆₀₀ of 5 U/ml. This takes ~20 min.
Total volume: _____ (_____ ml/pellet)
8. Once cells are completely resuspended, incubate for 45 min. on shaker in cold room, setting 4.
9. Pellet cells (15,000 x g, 10 min, 4°C).
10. Remove and **save the periplasmic fraction** (supernatant).
11. Remove 30 μ l of periplasmic fraction for SDS-PAGE analysis. Add 30 μ l sample buffer – freeze.
12. Resuspend the pellet in PBS (1/10 volume used during osmotic shock). Volume: _____.
13. Remove 10 μ l of the resuspended pellet for SDS-PAGE analysis. Add 10 μ l sample buffer – freeze.
14. Wash (5 min on rotator) 4 ml Ni beads with 40 ml pre-elution 10 buffer (300 mM NaCl, 50 mM Tris-Cl, pH 8.0, 10 mM imidazole, 10% glycerol), then pellet (10 min, 1,400 g, RT) and discard S/N wash.
15. Incubate in 250 ml conical bottles, the periplasmic fraction with the pre-washed Ni beads (2 hours, 4°C, on rotator). Parafilm the lids.
16. Pellet the Ni beads (10 min, 2,500 rpm, RT), remove and save the Unbound fraction. Combine 20 μ l of Unbound with 20 μ l sample buffer – freeze.
17. Add the Ni beads to a column and allow the beads to settle for 30 min.
18. Wash the beads with pre-elution buffer 10, until baseline is flat. Collect and analyze this wash by SDS-PAGE.
19. Elute the bound protein with elution buffer (300 mM NaCl, 50 mM Tris-Cl, pH 8.0, 250 mM imidazole, 10% glycerol). Collect the eluted protein in smallest volume possible – watch chart recording to ascertain when the majority of the protein has eluted from the column.
20. Take A₂₈₀ measurements of the eluted protein. A₂₈₀: _____ (= _____ mg/ml) blank = _____
Eluate volume: _____, Yield from Ni column: _____ mg.
21. Remove 10 μ l of eluted protein for SDS-PAGE analysis. Add 10 μ l sample buffer – freeze.
22. Remove and save the Ni beads, recharge as described below.
23. Concentrate the eluate using the AMICON spin cell concentrator and Millipore membrane #13612. Starting volume _____, ending volume _____.
24. A₂₈₀ post concentration _____ (1:100 with water). (= _____ mg/ml)
25. Filter through pre-rinsed (with 50 mM Tris, pH 8.0, 150 mM NaCl) 0.2 μ m filter.
26. A₂₈₀ post filtration _____ (1:100 with water). (= _____ mg/ml)
27. Load onto pre-equilibrated Superdex 75. Buffer used to wash and elute the Superdex 75 column: _____.
28. Analyze FPLC fractions by SDS-PAGE and pool the fractions of interest.
29. Concentrate the pooled fractions to desired concentration using either an AMICON spin cell concentrator/Millipore membrane or Amicon Ultra-4 centrifugal filter device (10,000 MWCOs).

Recharging Ni beads:

1. Wash Ni beads in 5 column volumes of nanopure water.
2. Wash Ni beads in 3 column volumes of strip buffer (100 mM EDTA, 500 mM NaCl, 20 mM Tris, pH 8.0). Beads should turn white.
3. Wash Ni beads in 10 column volumes of nanopure water.
4. Wash Ni beads in 5 column volumes of charge buffer (100 mM NiCl₂).
5. Wash Ni beads in 10 column volumes of nanopure water.
6. Wash Ni beads in 5 column volumes of pre-elution buffer 10. Store beads in pre-elution buffer 10 plus either 10% EtOH or NaN₃.